INFORMATION LEAFLET FOR DKF IMMUNOLOGICAL FECAL OCCULT BLOOD RAPID TEST

EN

DKF i(FOB) Rapid Test

DKF Immunological Fecal Occult Blood Rapid Test

INTENDED USES

The DKF i(FOB) Rapid Test is a rapid test for the qualitative detection of low levels of fecal occult blood in human feces, providing results within 5-10 minutes. It is intended for professional use as a screening test and aid in the detection of bleeding caused by a number of lower gastrointestinal disorders, e.g., diverticulitis, ulcerative colitis, colon polyps, colorectal carcinoma and Crohn's disease. This test provides a preliminary screening result; alternative test method(s) should be considered to confirm the test result obtained by this device.

SUMMARY

Many diseases can cause hidden blood in the feces which is known as fecal occult blood. In the early stages, lower gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood. The immunological Fecal Occult Blood i(FOB) Rapid Test is a rapid test for qualitative detection of low levels of fecal occult blood in the human feces within 5-10 minutes. The i(FOB) Rapid Test is designed to be less sensitive to dietary influences compared to older fecal occult blood tests. However, it is still important for patients to adhere to any dietary or medication guidelines provided by specialists prior to taking the test to ensure the most accurate results. Therefore, the test can be performed in approximately 10 minutes and does not require highly skilled personnel. It is important to note that a positive result from i(FOB) Rapid Test is an initial screening outcome and should be followed up with consultation from a healthcare professional for further investigation and diagnosis to determine the underlying cause and appropriate treatment.

TEST PRINCIPLE

The DKF i(FOB) Rapid Test is a lateral flow chromatographic immunoassay. The test contains special proteins known as antibodies that specifically recognize human hemoglobin, those antibodies are attached to tiny gold particles (gold-labeled antibodies), and placed on specific sites of the test strip. When you add the extracted stool sample to the sample well of the test device, if there's any hemoglobin present in the sample it will bind to the gold-labeled specific antibodies forming an immunocomplex. This immunocomplex then migrates along the test strip until it reaches the area precoated with the hemoglobin specific antibodies in which it will be caught and a pink line appears on T line of the cassette, indicating a positive result. If there is no hemoglobin in the stool, no pink line will appear in this area, indicating a negative result. There's also a control line C Line on the test that should always appear regardless of whether hemoglobin is detected or not. This line is important because it shows that the test has worked correctly. If this C control line doesn't appear, this means that the test hasn't worked properly, and you'll need to repeat the test using a new test device. This method is sensitive enough to detect very low levels of blood in the stool, providing a quick way to screen for conditions that might not show obvious symptoms initially.

MATERIALS NEEDED

- Test device individually foil pouched with a desiccant bag
- Extraction Tube with sample collection stick
- Solution Buffer vial (1.5 ml)
- Stool cup
- Clock or timer
- Clean and disposable gloves
- Instruction of Use (product leaflet)

SPECIMEN COLLECTION AND STORAGE CONDITION - FIGURE 1

This test is performed on fecal specimen. To collect this sample, proceed with the following protocol:

- 1- Collect sufficient quantity of fecal samples (1-2 g) in a clean, dry stool cup.
- 2- Open the extraction tube by unscrewing the cap, and transfer the solution buffer into the bottle.
- 3- Use the collection stick of the cap to collect about 50 mg of feces by piercing the stool sample in at least five different spots. Note: Do not scoop the fecal specimen. Ensure that stool specimen is only in the grooves of the collection stick. Excess stool specimen may lead to an invalid test result.
- 4- Place the collection stick into the extraction tube filled with buffer solution and then tighten the cap securely to close it.
- 5- Shake the extraction tube vigorously for 15-20 sec, to mix the specimen and the solution buffer to ensure a homogenous liquid suspension. Let the mixture sit for 2 minutes.

Specimen collected may be stored for 3 days at 2-8°C if not tested within 6 hours. For long term storage, specimens should be kept below -20°C. If not tested within 1 hour of preparation, specimens in the collection tube can be stored for up to 6 months at -20°C.

Specimen should not be collected during or within three days of a menstrual period, or if the patient is suffering from bleeding hemorrhoids, constipation-related bleeding, or urinary bleeding.

For laboratory personnel: Pack the specimens in compliance with applicable regulations for transportation of etiological agents, in case they need to be shipped.

Alcohol, aspirin and other medications taken in excess should be discontinued at least 48 hours prior to testing because it may cause gastrointestinal irritation resulting in occult bleeding.

Dietary restrictions are not necessary. Do not use the test if there is visible blood in the feces, as this is not 'occult' and requires direct medical assessment.

TEST PROCEDURE - FIGURE 2

If stored in a refrigerator or freezer, allow the specimen and test components to reach room temperature before use.


If the specimen was frozen, ensure it is thoroughly thawed and then mix well before performing the assay.

After preparing the specimen as previously described, promptly remove the test device from its sealed pouch and place it on a clean, flat surface for immediate use.

Hold the specimen collection tube upright, then snap off its tip. Invert the tube and gently squeeze to dispense 3 full drops of the extracted specimen into the specimen well (S) of the test device. Then, start the timer.

Avoid trapping air bubbles in the specimen well (S). Do not overload samples.

Wait for the appearance of pink line(s). Interpret the results between 5 and 10 minutes. Do not interpret the results after 10 minutes, as they may no longer be accurate.

Notes:

Only open the sealed pouch when you are ready to conduct the assay.

Best results will be obtained if the assays performed within 6 hours after collection. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

Results can be read at 5 to 10 minutes. Positive results may be visible as soon as 1 minute. However, to confirm negative results, wait until the full 10 minutes period has elapsed. However, Results interpreted after the 10 minutes mark are invalid and the test must be repeated for accurate results.

After interpreting the results, dispose of the used device in accordance with local regulations governing biomedical waste.

If the specimen does not migrate, likely due to the presence of particles, centrifuge the extracted specimen in the extraction buffer vial. Collect 80 µL of supernatant and dispense it into the specimen well (S) of a new test cassette and start following the instructions mentioned above.

RESULTS INTERPRETATION - FIGURE 3

Positive Result:

- Description: Two distinct pink lines appears, one in the control region (C) and another in the test region (T).
- Meaning: Indicates hHb (human hemoglobin) is present with detectable concentration in the tested sample, signifying a positive or reactive result.
- Note: The color intensity of the test line (T) varies based on the concentration of Fecal Occult Blood in the specimen. Any shade in the test region indicates a positive result.

Negative Result:

- Description: One pink line appears in the control region (C) with no colored line in the test region (T).
- Meaning: Indicates hHb isn't present or present with undetectable concentration in the tested sample, signifying a negative or non-reactive result.

Invalid Result:

- Description & Cause: Control line fails to appear, likely due to insufficient specimen volume or incorrect procedural techniques.
- Action: Review the procedure and repeat the test with a new device. If issues persist, stop using the kit and contact your local distributor. If excess feces are the issue, collect a new specimen and retest.

WARNINGS AND PRECAUTIONS

- Keep out of children's reach.
- Use before the expiration date and do not reuse.
- For in vitro diagnostic use only.
- Store and handle under appropriate temperature and humidity conditions.
- Do not use if the foil pouch is damaged or improperly sealed.
- Avoid substituting components with those from other kits.
- Avoid eating, drinking, or smoking where specimens or kits are handled.
- Follow standard biosafety guidelines and wear protective clothing.
- Treat all specimens as potentially infectious.
- Dispose of specimens and used devices in bio-hazard containers according to local, national, or regional regulations.
- Thoroughly read the package insert before use to avoid inaccurate results.

STORAGE AND STABILITY

- Avoid direct sunlight, moisture, and heat.
- Store or refrigerate between 2-30°C. Do not freeze.
- Use the test device before the expiration date and within one hour after opening the pouch.
- Ensure kit components are contamination-free. Do not use if contamination or precipitation is evident.

QUALITY CONTROL

Internal Control: Presence of the C line indicates proper test execution. If the C line doesn't appear, repeat the test with a new device.

External Control: Use external controls (positive & negative) under the following specific circumstances to ensure assay accuracy:

- A new lot of test kits is used.
- A new shipment of test kits is used.
- To investigate the cause of repeated invalid results.
- The temperature of the test area falls outside of 15-30°C.
- The temperature during storage of the kit falls outside of 2-30°C.
- To verify a higher-than-expected frequency of positive or negative results.
- A new operator uses the test kit, prior to performing testing of specimens.

TEST LIMITATIONS

- Inaccurate results may arise from failure to adhere to test procedures, including improper sample collection, handling, or storage.
- Humidity and extreme temperature conditions can adversely affect test results.
- This test is designated for in vitro diagnostic use only.
- This test detects human hemoglobin in fecal specimens. However, it should be noted that fecal blood may originate from sources other than colorectal bleeding.
- Like all diagnostic tests, results from this test must be interpreted in the context of other clinical information available to the physician. It is an aid in diagnosis and should not replace other diagnostic procedures, such as gastrointestinal (G.I.) endoscopy, colonoscopy, or radiographic evaluations.
- Other clinically available tests are required if questionable results are obtained.
- A positive test result should prompt further diagnostic investigation to determine the exact cause of occult blood in the feces.
- Negative results do not rule out gastrointestinal disorders, including polyps and colorectal cancer, as these conditions may not always lead to detectable levels of blood in the stool or may bleed intermittently.
- There is always a possibility that false results will occur due to the presence of interfering substances in the specimen or factors beyond the control of the manufacturer, such as technical or procedural errors associated with the testing.
- While the test is highly accurate in detecting fecal hemoglobin, false results are possible. Confirmatory testing is advised in case of questionable results. A definitive clinical diagnosis should not be based solely on this test but should consider all clinical and laboratory findings.

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