INFORMATION LEAFLET FOR DKF HELICOBACTER PYLORI STOOL ANTIGEN RAPID TEST

EN

DKF H.Pylori Stool Ag Rapid Test Helicobacter Pylori Stool Antigen Rapid Test

INTENDED USES

The DKF H. pylori Stool Ag Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of Helicobacter pylori antigens in human stool specimen. It is for in vitro diagnostic use only, and provides a preliminary test result to aid in the diagnosis of active H. pylori infection.

SUMMARY

Helicobacter pylori (H. pylori) is a gram-negative spiral-shaped bacterium that colonizes the gastric mucosa of up to 50% of the population worldwide, with a higher prevalence in developing countries. H. pylori infection is a risk factor for a variety of gastrointestinal diseases including non-ulcer dyspepsia, duodenal and gastric ulcers and active chronic gastritis. Therefore, the hypothesis that elimination of H. pylori may reduce the incidence of gastric cancer is a current area of research. H. pylori can be transmitted through the ingestion of food or water that is polluted with fecal matter. Antibiotics in combination with bismuth compounds have been shown to be effective in treating active H. pylori infection. H. pylori infection is currently detected by invasive testing methods based on endoscopy and biopsy (i.e. histology, culture) or non-invasive testing methods such as Urea Breath Test (UBT), serologic antibody rapid test and stool antigen rapid test. The stool antigen rapid test detects H. pylori antigens present in the feces, which indicates an active H. pylori infection compared to serologic antibody rapid tests which cannot distinguish between current active and past infections. Compared to UBT, H. pylori stool antigen rapid test doesn't require expensive lab equipment and use of a radioactive reagent. H. pylori Stool Ag Rapid Test can also be used to monitor the effectiveness of treatment and the recurrence of an infection, and it is not affected by the use of Proton Pump Inhibitors (PPI). The test can be performed within 10 minutes by minimally skilled personnel without the use of laboratory equipment.

TEST PRINCIPLE

The DKF H. pylori Stool Ag Rapid Test is a lateral flow chromatographic immunoassay. The test contains special proteins known as antibodies that specifically recognize H.pylori antigens, those antibodies are attached to tiny gold particles (gold-labeled antibodies), and placed on specific sites of the test strip. When you add a stool sample to the sample well of the test device, if there's any H.pylori antigens present in the sample it will bind to the gold-labeled specific antibodies forming an immunocomplex. This immunocomplex then migrates along the test strip until it reaches the area precoated with the H.pylori antigens specific antibodies in which it will be caught and a pink line appears on T line of the cassette, indicating a positive result. If there is no H.pylori antigens in the stool, no pink line will appear in this area, indicating a negative result. There's also a control line C Line on the test that should always appear regardless of whether H. pylori antigens are detected or not. This line is important because it shows that the test has worked correctly. If this C control line doesn't appear, this means that the test hasn't worked properly, and you'll need to repeat the test using a new test device. This method is sensitive enough to detect the active H.pylori infection.

MATERIALS NEEDED

- Test device individually foil pouched with a desiccant bag
- Extraction tube with sample collection stick
- Solution Buffer vial (2 ml)
- Product leaflet (user guide)
- Stool cup
- Clock or timer
- Clean and disposable gloves
- Plastic dropper (for watery stool only)

SPECIMEN COLLECTION AND STORAGE CONDITION - FIGURE 1 & 2

This test is performed on fecal specimen. To collect the fecal sample, proceed with the following protocol, If the fecal specimen is solid follow Procedure A, if the fecal specimen is watery follow Procedure B.

Procedure A: Solid fecal specimens (FIGURE 1)

- 1- Collect sufficient quantity of fecal samples (1-2 g) in a clean, dry stool cup.
- 2- Open the extraction tube by unscrewing the cap, and transfer the solution buffer into the bottle.
- 3- Use the collection stick of the cap to collect about 50 mg of feces by piercing the stool sample in at least five different spots.

Do not scoop the fecal specimen. Ensure that stool specimen is only in the grooves of the collection stick. Excess stool specimen may lead to an invalid test result.

- 4- Place the collection stick into the extraction tube filled with buffer solution and then tighten the cap securely to close it.
- 5- Shake the extraction tube vigorously for 15-20 sec, to mix the specimen and the solution buffer to ensure a homogenous liquid suspension then let the mixture sit for 2 minutes.

Procedure B: Watery fecal specimens (FIGURE 2)

- 1- Collect sufficient quantity of fecal samples (1-2 ml) in a clean, dry stool cup.
- 2- Open the extraction tube by unscrewing the cap, and transfer the solution buffer into the bottle.
- 3- Use the plastic dropper to transfer 2-3 drops (approximately 85-100 µL) of the watery fecal specimen from the stool cup to the extraction tube filled with buffer solution.
- 4- Replace the cap and tighten securely to close the tube.
- 5- Shake the extraction tube vigorously for 15-20 sec, to mix the specimen and the solution buffer to ensure a homogenous liquid suspension. Let the mixture sit for 2 minutes.

Specimen collected may be stored for 3 days at 2-8°C if not tested within 6 hours. For long term storage, specimens should be kept below -20°C. If not tested within 1 hour of preparation, specimens in the collection tube can be stored for up to 6 months at -20°C.

Specimen should not be collected during or within three days of a menstrual period, or if the patient is suffering from bleeding hemorrhoids, constipation-related bleeding, or urinary bleeding.

For laboratory personnel: Pack the specimens in compliance with applicable regulations for transportation of etiological agents, in case they need to be shipped.

Alcohol, aspirin and other medications taken in excess should be discontinued at least 48 hours prior to testing because it may cause gastrointestinal irritation resulting in occult bleeding.

Dietary restrictions are not necessary. Do not use the test if there is visible blood in the feces, as this is not 'occult' and requires direct medical assessment.

TEST PROCEDURE

Step 1: Bring the specimen and test components to room temperature if refrigerated. If frozen wait until the specimen is thawed, mix well prior performing the assay.

Step 2: When ready to test, open the pouch at the notch and remove the test device and place it on a clean, flat surface.


Step 3: Shake the stool extraction tube vigorously to ensure a homogenous liquid suspension.

Step 4: Hold the stool extraction tube vertically. Twist off the cap. Dispense 2-3 drops of the solution into the sample well of the cassette.

Do not overload the solution.

Step 5: Set the timer on 10 minutes.

Step 6: Results can be read at 10 minutes. Positive results can be visible in as short as 1 minute. Negative results must be confirmed at the end of the 15 minutes only. However, any results interpreted outside of the 10-15-minute window should be considered as invalid and must be repeated. Discard used device after interpreting the result following local requirements governing the disposal of device.

RESULTS INTERPRETATION

Positive Result: The appearance of two distinct pink lines indicates a positive result. One line should be visible in the control line region (C), and another pink line should be visible in the test line region (T).

Negative Result: A negative result is indicated by the presence of only one pink line in the control line region (C), with no line appearing in the test line region (T).

Invalid Result: If no colored line appears in the control line region (C), the test result is invalid, regardless of any line in the test line region (T). An invalid result suggests a possible issue with the test procedure or that the test device may be faulty. In such cases, it is recommended to repeat the test with a new device.

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic Use. Do not use after the expiration date & Read the package insert before performing the test.
- Bring all reagents to room temperature (15-30°C) before use.
- Do not open the sealed pouch, unless ready to perform the test.
- Avoid scooping out too much stool sample, as excessive fecal matter can clog the sample pad, interfering with sample migration.
- Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
- Do not smoke, drink or eat in areas where specimens or kit reagents are being handled.
- Avoid buffer solution contact with skin or eyes. Do not ingest.
- Dispose of all specimens and materials used to perform the test in accordance with bio-hazardous waste regulations.
- Follow local, state, and federal guidelines for the safe disposal of such materials. This typically involves placing used test components and specimens in a designated biohazard container.
- Consult your facility's waste management policies or local health authorities for specific disposal instructions relevant to your region.
- Consider any materials of human origin as infectious and handle them using standard biosafety procedures.
- The test results should be read 10-15 minutes after a specimen is applied to the sample well of the device. Any results interpreted outside of the 10-15 minutes window should be considered invalid and must be repeated.
- Do not perform the test in a room with strong air flow such as electric fan or strong air conditioning.

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unopened test devices at 2-30°C. If stored at 2-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable before it's expiration date. Do not freeze the kit or expose the kit to temperatures above 30°C.

QUALITY CONTROL

Internal Control: This test includes an internal control feature, the control line C line that should appear after the sample is added, indicating a correct functioning of the test. If the C line does not appear, review the entire procedure and repeat the test with a new device.


External Control: Control standards are not supplied with this kit. However, it is strongly recommended to test both positive and negative controls as part of good laboratory practice. This helps confirm the test procedure's accuracy and verify the proper performance of the test. External controls are particularly important in the following situations:

- When a new operator is using the kit for the first time.
- When starting a new lot of test kits.
- When opening a new shipment of test kits.
- If the storage temperature of the kits falls outside the range of 2-30°C.
- If the temperature of the testing area is outside the range of 15-30°C.
- To confirm the accuracy of an unusual high frequency of either positive or negative results.
- To investigate the cause of repeated invalid results.

TEST LIMITATIONS

- Strict Adherence to Procedure: Accurate results depend on strict adherence to the assay procedure and the results interpretation sections, particularly the specimen collection and storage condition. Deviations from the prescribed procedure may lead to inaccurate results.
- Qualitative Detection: The H. pylori Ag Rapid Test is designed for the qualitative detection of H. pylori antigen in human fecal specimens. The intensity of the test line does not necessarily correlate with the antigen concentration in the specimen.
- Negative Results and Infection Possibility: A negative or non-reactive result indicates the absence of detectable H. pylori antigen. However, a negative result does not rule out the possibility of H. pylori infection.
- Detection Limits and Non-Reactive Results: A negative or non-reactive result may occur if the level of H. pylori antigen in the specimen is below the detection limit of the assay or if the antigens are not present in the collected fecal specimen.
- Relation with Fecal Occult Blood (FOB) Tests: It is reported that approximately 39.3% of specimens with positive FOB test results also test positive for H. pylori. Thus, specimens testing positive with an FOB test may also yield a positive result with this rapid test.
- Further Testing if Symptoms Persist: If symptoms continue and the result from this rapid test is negative or non-reactive, alternative testing methods are recommended for further evaluation.
- Interpretation of Results: The results from this rapid test should be interpreted in conjunction with other diagnostic procedures and clinical findings to make a comprehensive assessment.

FIGURE 1

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