INFORMATION LEAFLET FOR DKF HEPATITIS C VIRUS ANTIBODY RAPID TEST

EN

DKF HCV Ab Rapid Test DKF Hepatitis C Virus Antibody Rapid Test

INTENDED USES

The Hepatitis C Virus Antibody Rapid Test is a single-use, rapid chromatographic immunoassay for the qualitative detection of antibodies against Hepatitis C Virus in human venous and capillary whole blood, serum, or plasma. This test is a screening tool that provides a preliminary result that must be confirmed using another diagnostic test. It is intended for in vitro diagnostic use only.

SUMMARY

HCV is the causative agent of hepatitis C, a liver infection that can be either acute or chronic. The virus is primarily transmitted through unsafe injection practices, inadequate sterilization of medical equipments, and the transfusion of unscreened blood and blood products. Individuals infected with HCV produce antibodies against this virus, and the presence of these antibodies in the blood indicates a current or past infection with HCV. This rapid test detects these antibodies using human serum, plasma, or whole blood samples. The test provides an initial screening result; further diagnostic tests should be performed to confirm the obtained result.

TEST PRINCIPLE

The DKF HCV Ab Rapid Test detects HCV antibodies in the sample. When a sample (whole blood, serum, plasma) is mixed with the provided buffer solution and applied to the sample well, it reacts with recombinant HCV antigen conjugated to a colored particle (colloidal gold) precoated on the sample pad of the test. This mixture migrates by capillary action to react with membrane-fixed reagents on the Test region (T) and then the Control region (C). If sufficient anti-HCV antibodies are present in the specimen, colored band will form at the Test region (T) of the membrane, indicating a positive result. If this band is absent, it indicates a negative result. A colored line should always appear at the Control line region (C) of the membrane, serving as a procedural control. It confirms that a proper volume of the specimen has been added, membrane wicking has occurred, and the Colloidal Gold conjugate is functional. The absence of the control band indicates an invalid test.

MATERIALS

- | | |
|---|---|
| • PROVIDED <ul style="list-style-type: none">. Test device individually foil pouched with a desiccant bag. Buffer Solution. Safety Sterile Lancet (for fingerstick whole blood). Capillary tube with mark line (for fingerstick whole blood). Product leaflet (user guide) | • NOT PROVIDED <ul style="list-style-type: none">. Alcohol pad or water and soap. Timer. Clean and disposable gloves. Centrifuge machine (for plasma and serum). Specimen collection container, dropper or pipette (for venous whole blood, serum or plasma) |
|---|---|

SPECIMEN COLLECTION AND STORAGE CONDITIONS - FIGURE 1

This HCV Ab Rapid Test can be performed using whole blood (from either venipuncture or fingerstick), plasma, or serum.

Follow the procedures below based on the specimen type:

• Whole Blood (Fingerstick) - FIGURE 1:

1. Wash the patient's hands with soap and warm water or clean the selected finger thoroughly with an alcohol pad. Allow it to air dry for 30 seconds.
 2. To increase blood flow, massage the hand from the wrist to the selected fingertip (preferably the middle or ring finger) without touching the puncture site.
 3. Use the provided safety sterile lancet to puncture the skin following the steps in the figure bellow. Wipe away the first drop of blood to remove any contamination.
 4. Encourage the formation of a rounded drop of blood by gently massaging the hand from wrist to palm to fingertip.
 5. Collect the blood using a capillary tube by touching its end to the blood drop. Ensure the tube fills to the marked line while avoiding air bubbles.
 - **Venous Whole Blood:** Collect the specimen into a collection tube containing an appropriate anticoagulant following standard venipuncture techniques.
 - **Serum:** Collect whole blood into a tube without anticoagulants. Centrifuge at 3000 rpm for at least 5 minutes to separate the serum. Use the supernatant for testing.
 - **Plasma:** Collect whole blood into a tube containing an anticoagulant. Gently invert the tube several times, then centrifuge at 3000 rpm for at least 5 minutes. Use the supernatant plasma for testing.
- For both serum and plasma,** separate them from the blood as soon as possible to prevent hemolysis.

Conditions:

- Use only clear, non-hemolyzed specimens for testing.
- Avoid leaving specimens at room temperature for extended periods.
- Perform tests immediately after collection when possible.
- Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For longer storage, keep specimens below -20°C.
- Before testing, bring specimens to room temperature. Thaw frozen specimens completely and mix well. Avoid repeated freezing and thawing (more than 3 times).
- Do not freeze whole blood specimens. Fingerstick whole blood should be tested immediately.
- If specimens need to be shipped, pack them in accordance with federal/local regulations for the transport of etiologic agents.

TEST PROCEDURE - FIGURE 2

1. Collect the sample as previously described.
 2. Put the test kit and sample at normal room temperature (15-30°C).
 3. Open the test pouch and take out the test cassette.
 4. Put the cassette on a clean, dry and flat surface.
 5. Add the sample to the sample well of the test cassette based on its type as follows:
For Whole Blood Sample from a Finger Prick:
 - After filling the capillary tube with blood until the mark, put it into the test cassette's sample area and press on its top to dispense the blood.For Serum, Plasma, and Venous Whole Blood Samples:
 - Hold the dropper straight up and put 1 drop of venous whole blood or 2 drops of serum /plasma into the cassette's sample well and press on its top to dispense the sample, or use a pipette as shown in the figure 2.
 6. Add 3 drops of buffer solution.
 7. Start the timer for 10 - 15 minutes.
 8. Wait for colored lines to show up and interpret the results. Don't read the results if it's been over 20 minutes since they might be incorrect.
- Notes:**
- Do not allow the specimen or buffer to enter the result viewing area.
 - It is critical to use the correct volume of buffer solution to ensure the test's validity.
 - If no migration is observed in the test window within one minute, add an additional drop of buffer to the sample well.
 - Avoid using the buffer if it has been more than 6 months since first opening the vial.
 - Positive results might appear as quickly as one minute for samples with a high viral load, but it is essential not to interpret the results before 10 minutes or after 20 minutes.

RESULTS INTERPRETATION - FIGURE 3


POSITIVE: Two colored bands appear on the membrane. One band appears in the control region (C), and another band appears in the test region (T). A positive result indicates the detection of HCV antibodies in the specimen.

Note: The intensity of the color in the Test line region (T) will vary depending on the concentration of HCV antibodies present in the specimen. Therefore, any shade of color in the test region should be considered positive.

NEGATIVE: Only one colored band appears, in the control region (C). No apparent colored band appears in the test region (T). A negative result indicates that HCV antibodies are not detected in the specimen.

INVALID: No band appears in the control region (C). Results from any test that does not produce a control band at the specified read time must be discarded. Possible causes of an invalid result include insufficient specimen volume or incorrect procedural techniques. If this occurs, review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit and contact your local distributor.

Note: It is normal for the control band to appear slightly lighter with very strong positive samples, as long as it is distinctly visible.

WARNINGS AND PRECAUTIONS


- Intended for in vitro diagnostic use only.
- Thoroughly read this package insert before performing the test. Failure to follow these instructions may give inaccurate test results.
- Do not use expired devices or buffers.
- Each device is designed for single use. Do not reuse.
- Do not mix or interchange components between different lots or different products.
- Open the sealed pouch only when ready to perform the assay. Do not use the test if the pouch is opened or damaged.
- Environmental conditions such as humidity and temperature can affect results. Perform the test in the prescribed conditions (15-30°C). Once opened, use the device promptly.
- Avoid performing the test in areas with strong airflow (e.g., near an electric fan or under strong air-conditioning).
- Do not eat, drink, or smoke in areas where specimens or test kits are handled.
- Avoid direct contact with the specimen pad, membrane, or result window of the test device to prevent inaccurate results.
- Dispose of all used specimens and testing materials as bio-hazardous waste.
- Treat all specimens as potentially infectious. Follow standard precautions against microbiological hazards throughout the procedure and adhere to proper disposal methods.
- Wear protective clothing and disposable gloves while handling kit reagents and clinical specimens. Wash hands thoroughly after testing.

STORAGE AND STABILITY

- Store the kit at room temperature or refrigerated (2-30°C). DO NOT FREEZE.
- The test cassette is stable until its expiration date.
- Keep the cassette in the sealed pouch and in a dry condition until use.
- Do not use the test after the expiration date or if the pouch is damaged.
- Do not reuse tests.
- Protect kit components from contamination. Avoid using the kit if there is evidence of microbial contamination or precipitation.
- Biological contamination of dispensing equipment or reagents can lead to false results.

QUALITY CONTROL

- Each test includes an internal procedural control: The appearance of a colored line in the control region (C) indicates a valid internal positive procedural control. It confirms sufficient specimen volume and correct procedural technique.
- External controls (positive and negative controls) are not included in this kit. However, testing positive and negative specimen controls is recommended as part of good laboratory practice. This ensures correct specimen volume, proper test procedure, and verifies the test device's performance.
- Handle the external negative and positive controls in the same manner as patient specimens.

TEST LIMITATIONS

Screening Use Only: This test is intended for screening purposes only. Positive results should be confirmed with alternative testing methods and clinical findings.

Indication of Antibodies Only: These tests indicate the presence of antibodies to HCV in serum, plasma, or whole blood specimens but should not be used as the sole criteria for diagnosing Hepatitis C viral infection.

Negative Results and Continued Symptoms: If test results are negative (non-reactive) but clinical symptoms persist, additional follow-up testing using other clinical methods is recommended. A negative result does not rule out the possibility of HCV infection, as it might occur if antibody levels are below the detection limit of the test or if antibodies are not present during the stage of disease when the sample is collected.

Interference Factors: In multiple blood transfusion recipients, there is a risk of weak reactive results in HCV negative samples.

Specimen Requirements: The hematocrit level of whole blood samples should be between 25% and 65%. Only clear, fresh, fluid whole blood, serum and plasma samples should be used. Fresh samples are preferred, but frozen thawed samples (except whole blood) can be used. Avoid shaking the sample.

In Vitro Diagnostic Use: These tests are for in vitro diagnostic use only and cannot determine the quantitative value or the rate of increase in HCV antibodies.

Note: All results must be interpreted in conjunction with other clinical information by a healthcare professional.

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DIAGNOSTIC KITS FACTORY S.A.R.L.