INFORMATION LEAFLET FOR DKF HUMAN IMMUNODEFICIENCY VIRUS RAPID TEST

EN

DKF HIV Ab Rapid Test

DKF Human Immunodeficiency Virus Antibody Rapid Test

INTENDED USES

The DKF HIV Ab Rapid Test is a quick, user-friendly test designed to detect the presence of antibodies against the Human Immunodeficiency Virus (HIV) types 1 and 2 in human serum, plasma, or whole blood from a finger prick or when drawn with anticoagulants. This test is intended for in vitro (outside the body, in a controlled environment) use to help in the diagnosis of HIV infection, a critical step in the management and prevention of AIDS (Acquired Immuno-Deficiency Syndrome).

SUMMARY

HIV is the virus that causes AIDS, a condition that weakens the immune system, making it hard for the body to fight off infections and diseases. There are two main types of HIV: HIV-1, which is more prevalent and the primary cause of HIV infections worldwide, and HIV-2, which is less common and mostly found in West Africa. Both types spread through contact with infected blood, sexual contact, and from mother to child during pregnancy, childbirth, or breastfeeding. Early and accurate detection of HIV is crucial for effective management and treatment. The HIV 1/2 Rapid Test Cassette offers a convenient method to quickly check HIV antibodies in the blood, serum, or plasma. These antibodies are the body's response to an HIV infection. By using a combination of proteins, that can specifically identify the presence of antibodies to both HIV-1 and HIV-2. Understanding the outcome of this test is essential. A positive result indicates that HIV antibodies were detected, and it's important to follow up with a healthcare provider for confirmatory testing and guidance. A negative result means that HIV antibodies were not detected in the sample. However, if recent exposure is suspected, re-testing after the window period (the initial few weeks post-exposure when antibodies may not be detectable) is recommended to ensure accuracy. This rapid test is a vital tool in the early detection of HIV, offering individuals a quick way to take charge of their health. However, it's important to remember that it should be used as part of a comprehensive approach to HIV testing and management, including confirmatory testing and professional medical advice.

TEST PRINCIPLE

This HIV 1/2 Ab Rapid Test is a qualitative, membrane-based immunoassay designed for detecting HIV antibodies in serum, plasma, or whole blood samples. This test operates on the principle of antigen-antibody reaction and chromatographic migration. The test cassette has a nitrocellulose membrane strip pre-coated with recombinant HIV antigens in the test line region T line. When a sample (blood, serum, or plasma) is applied to the sample pad, it first interacts with specific particles present on the membrane. This sample-particle mixture then moves along the membrane by capillary action towards the other end. As this mixture migrates, it encounters the recombinant HIV antigens. If antibodies to HIV-1 or HIV-2 are present in the sample, they bind to these antigens, resulting in the appearance of a colored line in the test line region, indicating a positive result. Conversely, if no HIV antibodies present in the sample, no colored line will appear in the test line region, indicating a negative result. To ensure the test has been performed correctly, a separate colored line will appear in the control region C line, confirming that the test is valid. This control line indicates that an adequate volume of the specimen has been added and that the specimen has migrated properly across the membrane. This simple yet effective mechanism allows for the rapid detection of HIV antibodies, facilitating diagnosis and management of HIV infection.

MATERIALS

NOT PROVIDED

- . Alcohol pad or water and soap
- . Timer
- . Clean and disposable gloves
- . Centrifuge machine (for plasma and serum)
- . Specimen collection container, dropper or pipette (for venous whole blood, serum or plasma)

PROVIDED

- . Test device individually foil pouched with a desiccant bag
- . Buffer Solution
- . Safety Sterile Lancet (for fingerstick whole blood)
- . Capillary tube with mark line (for fingerstick whole blood)
- . Product leaflet (user guide)

SPECIMEN COLLECTION AND STORAGE CONDITIONS - FIGURE 1

This HIV Ab rapid test can be performed using whole blood (from either venipuncture or fingerstick), plasma or serum. Follow the procedures below based on the specimen type:

• Whole Blood (Fingerstick) - FIGURE 1:

1. Wash the patient's hands with soap and warm water or clean the selected finger thoroughly with an alcohol pad. Allow it to air dry for 30 seconds.
2. To increase blood flow, massage the hand from the wrist to the selected fingertip (preferably the middle or ring finger) without touching the puncture site.
3. Use the provided safety sterile lancet to puncture the skin following the steps in the figure bellow. Wipe away the first drop of blood to remove any contamination.
4. Encourage the formation of a rounded drop of blood by gently massaging the hand from wrist to palm to fingertip.
5. Collect the blood using a capillary tube by touching its end to the blood drop. Ensure the tube fills to the marked line while avoiding air bubbles.

• **Venous Whole Blood:** Collect the specimen into a collection tube containing an appropriate anticoagulant following standard venipuncture techniques.

• **Serum:** Collect whole blood into a tube without anticoagulants. Centrifuge at 3000 rpm for at least 5 minutes to separate the serum. Use the supernatant for testing.

• **Plasma:** Collect whole blood into a tube containing an anticoagulant. Gently invert the tube several times, then centrifuge at 3000 rpm for at least 5 minutes. Use the supernatant plasma for testing.

For both serum and plasma, separate them from the blood as soon as possible to prevent hemolysis.

FIGURE 1

- Use only clear, non-hemolyzed specimens for testing.
- Avoid leaving specimens at room temperature for extended periods.
- Perform tests immediately after collection when possible.
- Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For longer storage, keep specimens below -20°C.
- Before testing, bring specimens to room temperature. Thaw frozen specimens completely and mix well. Avoid repeated freezing and thawing (more than 3 times).
- Do not freeze whole blood specimens. Fingerstick whole blood should be tested immediately.
- If specimens need to be shipped, pack them in accordance with federal/local regulations for the transport of etiologic agents.

TEST PROCEDURE - FIGURE 2

1. Collect the sample as previously described.
2. Put the test kit and sample at normal room temperature (15-30°C).
3. Open the test pouch and take out the test cassette.
4. Put the cassette on a clean, dry and flat surface.
5. Add the sample to the sample well of the test cassette based on its type as follows:
For Whole Blood Sample from a Finger Prick:
 - After filling the capillary tube with blood until the mark, put it into the test cassette's sample area and press on its top to dispense the blood.For Serum, Plasma, and Venous Whole Blood Samples:
 - Hold the dropper straight up and put 1 drop of venous whole blood or 2 drops of serum /plasma into the cassette's sample well and press on its top to dispense the sample, or use a pipette as shown in the figure 2.
6. Add 3 drops of buffer solution.
7. Start the timer for 10 - 15 minutes.
8. Wait for colored lines to show up and interpret the results. Don't read the results if it's been over 20 minutes since they might be incorrect.

Notes:

- Do not allow the specimen or buffer to enter the result viewing area.
- It is critical to use the correct volume of buffer solution to ensure the test's validity.
- If no migration is observed in the test window within one minute, add an additional drop of buffer to the sample well.
- Avoid using the buffer if it has been more than 6 months since first opening the vial.
- Positive results might appear as quickly as one minute for samples with a high viral load, but it is essential not to interpret the results before 10 minutes or after 20 minutes.


FIGURE 2

RESULTS INTERPRETATION - FIGURE 3

POSITIVE: Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T). A positive result in the Test region indicates the detection of HIV antibodies in the specimen. Note: The intensity of the color in the test line region (T) will vary depending on the concentration of HIV antibodies present in the specimen. Therefore, any shade of color in the test region should be considered positive.

NEGATIVE: Only one colored band appears in the control region (C). No apparent colored band appears in the test region (T). A negative result in the Test region indicates the absence of HIV antibodies in the specimen.

INVALID: The control band fails to appear. Results from any test that has not produced a control band at the specified read time must be discarded. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

Note: It is normal to observe a slightly lightened control band with very strong positive samples, as long as it is distinctly visible.

FIGURE 3

Positive+

Negative -

Invalid

STORAGE AND STABILITY

This Test is designed for storage at room temperature or in a refrigerated environment, spanning from 2°C to 30°C. Each test cassette remains stable and should perform accurately up until its expiration date. It is crucial that the test cassettes are kept sealed within their pouches and stored in dry conditions until they are ready to be used. Freezing the test kit is not recommended and can compromise the integrity and accuracy of the test. Test cassettes should not be used if they have passed their expiration date, or if the pouch has been damaged in any way. Reusing test cassettes is also not advised. Care must be taken to avoid contamination of the kit's components. The kit should not be used if there are signs of microbial contamination or if precipitation is observed, as these conditions can lead to false results. The integrity of dispensing equipment and reagents should be maintained to ensure accurate test outcomes.

WARNINGS AND PRECAUTIONS

- This test is for in vitro diagnostic use only.
- Carefully read and follow the instructions provided with the test kit. Failure to adhere to the instructions may result in inaccurate test results.
- Do not use the test devices or buffer solutions past their expiration dates. Ensure that all components are stored under recommended conditions and are within their usage period.
- Each test device is designed for a single use. Do not attempt to reuse test cassettes.
- Do not mix components from different kit lots or interchange components with those from other products, as this could compromise test results.
- Open the test device pouch only when ready to perform the test. Inspect the pouch for any damage before opening, and do not use the test if the pouch is damaged.
- Perform the test in an environment with a temperature range of 15-30°C to ensure accurate results. Avoid testing in areas with strong airflow, such as near fans or in heavily air-conditioned rooms.
- Do not eat, drink, or smoke in the area where the tests are being conducted. Ensure that the testing area is clean and well organized.
- Direct contact with the specimen pad, membrane, or result window should be avoided to prevent contamination and ensure accurate results.
- Treat all specimens as potentially infectious. Follow established biosafety protocols for handling and disposal of clinical specimens and testing materials.
- Test cassettes and any other used test materials should be disposed of according to local regulations for biohazardous materials. Autoclaving before disposal is recommended to ensure biological safety.
- Wear appropriate protective clothing and disposable gloves while performing the test and handling specimens. Wash hands thoroughly after the test is completed.

QUALITY CONTROL

Quality control measures are integral to this test to ensure the reliability of test results. An internal procedural control is incorporated within the test, manifesting as a colored line in the control region (C) on the test cassette. The appearance of this line is an indication of a valid test, confirming that an adequate volume of specimen has been applied and that the procedural technique has been correctly followed.

Although external controls are not included within this kit, it is highly recommended to use both positive and negative control specimens as part of routine laboratory practice. Testing control specimens helps to verify the accuracy of the test device by ensuring that the specimen volume is sufficient and that the procedural techniques used are correct. Handling these controls in the same manner as patient specimens is essential for maintaining the integrity of the testing process.

Implementing both internal and external quality control measures helps to maintain the accuracy and reliability of the rapid test, thereby ensuring that test results can be trusted for clinical decision-making.

TEST LIMITATIONS

- This test may not detect HIV antibodies during the "window period" of infection, which is the time between the initial infection and the appearance of detectable antibodies. This period can vary from person to person. Consequently, a negative result does not preclude the possibility of HIV infection. If recent exposure is suspected, re-testing at a later date or confirmatory testing using different methodologies might be necessary.
- Although the test is designed to be highly specific to HIV-1 and HIV-2 antibodies, cross-reactivity with antibodies from other diseases or conditions can occasionally occur, leading to false-positive results. Confirmatory testing using an alternative method is recommended in case of a positive result.
- This test may not detect all subtypes or variants of HIV. While the test is designed to identify a broad range of HIV-1 and HIV-2 antibodies, novel or rare strains might not be detected.
- Incorrect performance of the test or failure to follow procedural instructions precisely can lead to errors in interpretation. This includes inadequate specimen volume, improper handling or timing, and environmental factors not adhering to recommended conditions.
- Individuals with compromised immune systems, such as those with advanced HIV disease, have delayed or diminished production of antibodies, potentially leading to false-negative results.
- Certain substances or treatments, including antiretroviral therapy, might influence the accuracy of the test by affecting the immune response or the detection of antibodies.
- This test is a screening tool and not a definitive diagnostic test. Positive results require confirmation by more specific tests, such as Western Blot or PCR-based methods, which can differentiate between HIV-1 and HIV-2 and confirm the presence of viral RNA.
- This test is not approved for screening blood or plasma donors. More sensitive laboratory tests are available and should be used for screening donated blood products.

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