INFORMATION LEAFLET FOR DKF OVULATION (LH) RAPID TEST

EN

DKF Ovulation (LH) Rapid Test

DKF Ovulation Luteinizing Hormon Rapid Test

INTENDED USES

The DKF Ovulation LH Rapid Test is a rapid, visual chromatographic immunoassay. It is a one-step screening test designed for the qualitative detection of human luteinizing hormone (LH) in human urine specimens to predict the timing of ovulation. It is intended for self-testing and in vitro diagnostic use only.

SUMMARY

Luteinizing hormone (LH) plays a crucial role in the regulation of ovulation and ovarian functions during a woman's menstrual cycle. Ovulation, the release of an egg from the ovary, occurs in the middle of the menstrual cycle. For pregnancy to occur, fertilization by sperm must happen within 24 hours of the egg release. Immediately before ovulation, the body produces a significant amount of LH "LH surge", triggering the release of a matured egg. This test detects the LH surge, indicating potential ovulation within the next 24-36 hours, thus identifying peak fertility periods. Note that an LH surge and ovulation may not occur in all cycles. For an accurate and personalized interpretation, especially for individuals with irregular cycles or those finding it difficult to interpret their results, consultation with healthcare professionals is highly recommended.

TEST PRINCIPLE

The test operates on a rapid chromatographic immunoassay for the qualitative detection of LH in urine, aiding ovulation detection. It employs a combination of antibodies, to selectively identify elevated LH levels. A colored band forms at the test region if LH levels are sufficient, indicating a positive result. The absence of this band suggests a negative result. A separate colored line in the control region acts as a procedural control, confirming the test's proper execution.

MATERIALS

- Provided:
 - Test device individually foil pouched with dropper and desiccant bag
 - Insert package
- Not provided:
 - Specimen collection container
 - Timer

WHEN TO START TESTING

Determining the optimal time to start ovulation testing is crucial for accurate results. Follow these steps to identify the correct day to begin testing:

- 1.Determine the Menstrual Cycle Length: The menstrual cycle length is the duration from the first day of menstruation (the start of a new period) to the day before the next period begins.
- 2.Calculate Days to Start Testing: After establishing the menstrual cycle length, use the chart below to find out how many days after your period you should begin testing.
- Locate your menstrual cycle length in the first row of the chart.
- Find the corresponding number of days to wait before starting testing in the second row.
- 3.Identify the Starting Day for Testing: Count the days, starting from the first day of your last period, based on the number obtained in the previous step. This count includes the first day of your last period.
- 4.Testing Frequency: It's recommended to test once daily over a span of five days to adequately capture the LH surge. You can perform the test at any time of the day but remember to repeat the test nearly at the same time for five days. Try not to urinate for four hours before your ovulation test, and avoid drinking too much liquid before the test.

Menstrual cycle length	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	More than 40
Days to count ahead	6	6	7	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	Consult a doctor

Example:

•If your usual menstrual cycle is 28 days, the chart suggests to start testing 11 days after your period begins. If your last period started on the 3rd, begin counting 11 days from the 3rd, and start testing on the 13th then continue testing until the 17th.

Special consideration:

- If your menstrual cycle is shorter than 21 days or longer than 40 days, it's advisable to consult a doctor.
- If unsure about your cycle length, consider starting your tests 11 days after the onset of your period, assuming an average cycle length of 28 days. Continue testing daily until you detect the LH surge.
- If the length of your menstrual cycle is uncertain or consistently varies, use the shortest cycle (21 days) as a reference for the chart. This adjustment might necessitate extending the testing period beyond five days.


SPECIMEN COLLECTION AND STORAGE CONDITIONS

This rapid test is intended for use with human urine specimens only.

- Restrict fluid intake approximately 2 hours before collecting your urine to avoid diluting the sample.
- Determine the optimal time for urine collection. For best results, collect urine at approximately the same time each day, ideally between 10 AM and 8 PM, as LH synthesis occurs during these hours. Avoid using the first urine after waking up.
- Use a clean, dry container for urine collection.
- A fresh urine specimen is preferred. However, if testing is delayed, the specimen can be stored at room temperature for up to 8 hours or refrigerated at 2-8°C for up to 24 hours. Do not freeze. Allow refrigerated urine to reach room temperature before testing (about 30 minutes).
- If specimens need to be shipped, ensure they are packaged in compliance with all applicable transportation regulations for etiological agents.
- Avoid using turbid or contaminated specimens. Do not shake the container. Allow any sediment to settle and use only urine from the top of the container.
- Upon determining the start day for testing, collect your urine daily, recording the date, cycle day, and collection time. Refer to the TEST RESULTS CHART provided at the end of this package insert for recording these details.

TEST PROCEDURE

- 1.Bring the test and urine specimen to room temperature (15-30°C) before testing.
- 2.Open the pouch at room temperature, remove the test cassette, and use it promptly.
- 3.Place the test cassette on a clean, level surface. Using the dropper, transfer 3 full drops of urine into the specimen well (S) of the test cassette, avoiding air bubbles. Begin timing as per the instructions below.
- 4.Wait for colored lines to develop. Interpret the results at 5 minutes, but to confirm negative results, the full 10-minute reaction time is necessary. And do not read results after 10 minutes as they may not be reliable.

RESULTS INTERPRETATION

POSITIVE: Two distinct lines are visible, and the line in the test line region (T) is as dark as or darker than the line in the control line region (C). This result indicates a probable LH surge, suggesting that ovulation is likely to occur within the next 24-36 hours.

NEGATIVE: Two lines are visible, but the line in the test line region (T) is lighter than the line in the control line region (C), or there is no line at all in the test line region (T). This result indicates that no LH surge has been detected at this time.

INVALID: The control line does not appear. This outcome usually results from insufficient specimen volume or incorrect procedural techniques. Review the procedure and repeat the test using a new device. If the problem persists, stop using the test kit and contact your local distributor for assistance.

NOTES:

1. Variability in Test Line Intensity: The intensity of the color in the test line region (T) can vary, reflecting the concentration of LH in the specimen. A faint test line (T), indicating a negative result, means that LH is present but not at the surge level.
2. Variation in Ovulation Timing: Not all women ovulate in the middle of their menstrual cycle. Therefore, it's possible not to observe a positive result during the chosen testing days for some users.
3. Common Causes of Invalid Results: The most likely reasons for failure to observe the control band include insufficient specimen volume, incorrect procedural execution, or the use of an expired test.

WARNINGS AND PRECAUTIONS

- This product is intended for professional in vitro diagnostic use only.
- Ensure you read and understand the entire testing procedure before initiating the test.
- Do not use the test kit past its expiration date. If the pouch is damaged or has been opened, do not use the test. Tests are designed for single use only and should not be reused.
- Avoid eating, drinking, or smoking in the area where specimens and test kits are handled to prevent contamination.
- The test kit must remain in its sealed pouch until it is ready to be used.
- Prevent cross-contamination by using a new, clean specimen collection container for each urine sample.
- Treat all specimens as potentially infectious. Follow standard precautions to avoid microbiological hazards throughout testing and adhere to proper disposal methods for specimens.
- Dispose of used tests in accordance with local regulations concerning biomedical waste.
- Be aware that environmental conditions such as humidity and temperature can impact test results. Store the test kit in a dry location at temperatures between 2-30°C. Do not freeze the test kit.

STORAGE AND STABILITY

- Store the test kits in their original sealed pouches at room temperature or refrigerated, within the range of 2-30°C (36-86°F).
- The test kits are stable until the expiration date. Do not use the test kits beyond this date.
- It is crucial to keep the test kit sealed in its pouch until ready for use to maintain its integrity.
- Freezing the test kit is strictly prohibited as it can compromise the effectiveness of the test.
- Protect the test kit from direct sunlight, moisture, and excessive heat to ensure accurate results.

QUALITY CONTROL

A Quality control checks are built into the testing process to ensure accuracy and reliability. A control line (C) in the control line region serves as an internal procedural control. The appearance of a colored line in this region verifies that sufficient specimen volume has been collected, an adequate membrane wicking has occurred and the test procedure has been correctly followed.

TEST LIMITATIONS

1. Professional Use
- The test is designed for in vitro diagnostic use by healthcare professionals.
2. Not a Contraceptive Method
- This test predicts ovulation to aid in conception planning. It should not be used as a method of contraception.
- Results from this test cannot assist in any form of contraception planning.

3. Safety Precautions

- Due to small parts and the nature of the test, it should be kept out of reach of children to prevent accidental ingestion or misuse.
- 4. Usage Restrictions
 - The test is not designed to provide accurate results for individuals who are pregnant, have reached menopause, or are undergoing hormone treatment, as these conditions significantly alter hormone levels and can skew results.
 - While common drugs such as pain relievers and antibiotics generally do not affect the test, medications containing human chorionic gonadotropin (hCG) or LH can lead to inaccurate results. Therefore, individuals taking such medications should avoid using this test.
 - Those on birth control pills, pregnant women, and individuals in menopause should not use the test, as it will not provide meaningful results under these conditions.

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