INFORMATION LEAFLET FOR DKF MULTIDRUG RAPID TEST

EN

DKF Multi-drug Rapid Test

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INTENDED USES

The Multi-Drug Rapid Test is a rapid, qualitative lateral flow chromatographic immunoassay for the screening of any combination of the following drugs and/or their metabolites in urine at or above the specified cut-off levels, the test is designed for simplicity and can be interpreted visually without specialized equipment. The speed and sensitivity of immunoassay have made it the most widely accepted method to screen for the presence of multiple drugs of abuse in urine, the detected drugs with their Cut-off concentrations are listed in the table below (Cut-off Levels are based on international standards or guidelines):

Abbreviation	Test	Cut-off
AMP	Amphetamine	500 ng/ml
MET	Methamphetamine	500 ng/ml
MOP	Morphine	300 ng/ml

This test provides only a preliminary result. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

AMP (Amphetamine) - Amphetamine is a central nervous system stimulant that affects chemicals in the brain and nerves contributing to hyperactivity and impulse control. Although used in treating attention deficit hyperactivity disorder (ADHD), its detection in human urine is widely used to assess amphetamine abuse. Effects include alertness, wakefulness, increased energy, reduced hunger, and an overall feeling of well-being. Overdose and extended usage can lead to severe and/or permanent damage to the nervous system. Amphetamines can be detected in urine within 3 hours after administration and can remain detectable for about 24-48 hours after the last dose. MET (Methamphetamine) - Methamphetamine, derived from its parent drug amphetamine, is a more potent central nervous system stimulant with longer-lasting effects. It is medically indicated for ADHD treatment and as a short-term component of weight-loss treatments but is highly addictive. Overdosage can cause severe symptoms including restlessness, confusion, and potentially fatal conditions like circulatory collapse and coma. Up to 43% of a dose is eliminated as amphetamine. MET concentrations can be significantly higher in abusers. MET's long-lasting effects contribute to its high potential for addiction and abuse. MOP (Morphine) - Morphine, a widely used analgesic, is effective for moderate to severe pain management. It acts directly on the central nervous system (CNS) to induce analgesia. However, repeated administration can lead to dependence and tolerance. Adverse effects include respiratory depression, drowsiness, and, in cases of overdose, coma and death. Morphine is metabolized extensively, with only 2-12% excreted as unchanged morphine in the urine. It is also a common metabolite of opiates like codeine and heroin, making morphine a key marker for opiate abuse.

TEST PRINCIPLE

The Multi-Drug Rapid Test is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a urine specimen migrates upward across the membrane by capillary action. If a drug is present in the urine specimen below its cutoff concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will appear in the corresponding test line region. If a drug is present in the urine specimen below its cutoff concentration, it will not saturate the binding sites of its specific antibody coated on the particles, and the colored line will not form in the corresponding test region. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper specimen volume has been added, correct procedural techniques were employed and membrane wicking has occurred.

MATERIALS

- PROVIDED
- Test panels sealed in an aluminum foil pouch with a desiccant bag
- Product leaflet (user guide)
- Not PROVIDED
- Specimen collection container
- Timer

SPECIMEN COLLECTION AND STORAGE CONDITION

The urine specimen must be collected in a clean and dry container, which should be securely closed and sterile. Urine collected at any time of the day may be used. If urine specimens exhibit visible precipitants, they should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

For optimal results, it is recommended to test urine specimens immediately after collection. However, if immediate testing is not feasible, the following storage conditions can be applied:


- Urine specimens can be stored at 15-30°C for up to 8 hours without significant loss of testing integrity.
- For storage extending beyond 8 hours, specimens should be refrigerated at 2-8°C and tested within 3 days.
- For long-term storage, specimens may be frozen at -20°C or lower. It is important to thaw and mix frozen specimens thoroughly before testing.

Please note: These storage guidelines apply specifically to urine specimens. The test kit itself should not be frozen, as freezing can damage the reagents and compromise test performance.

TEST PROCEDURE - FIGURE 1

If refrigerated, allow urine specimen to equilibrate to room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test panel from the sealed pouch and use it as soon as possible.
2. Collect urine specimen
3. Remove the cap. With the arrow pointing toward the urine specimen, immerse the test panel vertically in the urine specimen for at least 10 to 15 seconds. Immerse the dipstick to at least the level of the wavy lines, but not above the arrow on the test panel.
4. Put the cap back on the test panel and place it on a non-absorbent flat surface.
5. Start the timer and wait for the colored line(s) to appear.
6. The drug strip result should be read at 5 minutes. Do not interpret the result after 10 minutes

RESULTS INTERPRETATION


Do not read test results after 10 minutes.

The T Line should always be interpreted independently of the C Line. Do not compare color intensity of one test to another.

NOTE: The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

Negative Result: Two distinct lines will appear. One line will be located in the control region (C), indicating the test is valid. The second line, appearing in the test region (T), confirms a negative result for the drug. The color intensity of the test line may vary, but any shade of pink or red, however faint, should be considered negative.

Positive Result: A single line appears in the control region (C), indicating the test is valid. The absence of a line in the test region (T) signifies a positive result for the drug being tested.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use.
- Do not use after expiration date.
- Do not reuse the same test panel.
- If the pouch is damaged, discard the test.
- Do not open the sealed pouch, unless ready to conduct the test.
- The instructions must be followed exactly to obtain accurate results.
- Ensure that the test cup is properly sealed following specimen collection.
- Dispose all specimens and used assay materials as potentially bio-hazardous.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- Avoid cross-contamination by using new gloves and a new specimen collection container for each test. Ensure that the testing area is clean and free from residues of previous tests.

STORAGE AND STABILITY

The test kit should be stored as packaged in the sealed pouch, either at room temperature or refrigerated (2-30°C). Do not freeze the test kit. Freezing can damage the reagents within the test panel, potentially leading to inaccurate results. The test panel is stable through the expiration date marked on the sealed pouch. It is crucial to keep the test panel in the sealed pouch until ready for use. Do not use the test kit beyond its expiration date.

Clarification: The warning against freezing is specific to the test kit and its components. Urine specimens, in contrast, may be frozen for long-term storage as detailed in the "SPECIMEN COLLECTION & STORAGE CONDITIONS" section, provided they are properly thawed and mixed before testing.

QUALITY CONTROL

An internal procedural control is built inside the test. A colored line appears in the control region (C) is considered an internal positive procedural control, and indicates a valid test result, it confirms sufficient specimen volume and correct procedural technique. External controls are not supplied with this kit; however, it is recommended to test positive and negative specimen controls as a good laboratory practice to confirm sufficient specimen volume and the right operation procedure and to verify proper test performance of the test device. Handle the negative and positive controls in the same manner as patient specimens. It is recommended to test positive and negative external controls with each new batch of tests or at least once a month to ensure the test kit is functioning correctly and to verify the test procedure's accuracy. This practice helps to maintain the integrity of test results.

TEST LIMITATIONS

1. This test is designed for testing human urine only. This test provides only a preliminary qualitative analytical result. A more specific alternate method must be used in order to obtain a confirmed result.
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants such as bleach and/or alum or other strong oxidizing agents, in urine specimens may produce erroneous results regardless of the test method used. When suspected, the test should be repeated with another fresh urine specimen.
4. A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.

Cross-Reactivity and False Results:

- Medications: Certain prescription and over-the-counter medications may interfere with test results. For instance, medications containing pseudoephedrine might cause false-positive results for methamphetamine.
- Dietary Substances: Natural substances in foods, such as poppy seeds, have been known to cause false-positive results for opiates due to their morphine content.
- Interpreting Results: It's essential to consider an individual's history of prescription medication use, dietary intake, and other potential interfering substances when interpreting test results.

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